

ProtectIR

Accurate identification of solid and liquid threats

Extensive Identification Capabilities


Confidently identify thousands of chemical threats anywhere with the portable ProtectIR. Compact and lightweight, it's the ideal "grab and go" platform for rapidly assessing chemical threats.

Identify 23,000+ Solids and Liquids:

- Hazardous Chemicals
 - TICs
 - TIMs
 - VOCs
- Narcotics, including fentanyl
- Explosives

Features:

- Clear, visible results in seconds
- Simple operation
- 24/7/365 Reachback
- Designed for use in harsh environments
- Built-in wireless communications



Specifications

Technology:	FTIR (Fourier Transform Infrared Spectroscopy)
Sample Type:	Solids, Liquids
Sample Interface:	Diamond ATR (Attenuated Total Reflectance)
Size:	6.5" × 9.25" × 3" (16.51cm × 23.5 cm × 7.62cm)
Weight:	5 lbs (2.27 kgs)
Frequency Range:	4000–650 cm
Resolution:	4 cm
Display:	4.3", 480 × 272
Connections:	USB 2.0 WiFi (802.11 G)
Durability:	IP67, Mil spec 810G
Search Technique:	Automated Mixture Search
Alarms:	Protein, Explosive, CWA, Water (as an interferent)
Libraries:	23,000+ Powders and Liquids. TICs, Narcotics, CWA, Explosives, Consumer Products, Chemicals
Power:	Rechargeable Li Ion battery, 4 hour run time Line power transformer, 100–240 VAC, 50/60 Hz
Database:	ThreatAssist; integrated chemical information database
Decon:	Spray, wipe-down or dunk with bleach or detergent
Warranty:	1 Year Warranty 24/7/365 Support Reachback with data analysis
Remote Data Access:	Team Leader App available for download on iOS, Android, and Windows devices